In the Specification:

Please amend the specification as shown:

Please delete the paragraph on page 156, lines 7-20, and replace it with the following paragraph:

EXAMPLE 47

GSK3-B/Aurora Kinase Inhibitory Activity Assay

AuroraA (Upstate Discovery) or GSK3-β (Upstate Discovery) are diluted to 10nM and 7.5nM respectively in 25mM MOPS, pH 7.00, 25mg/ml BSA, 0.0025% Brij-35, 1.25% glycerol, 0.5mM EDTA, 25mM MgCl₂, 0.025% β-mercaptoethanol, 37.5mM ATP and and 10 μl mixed with 10 μl of substrate mix. The substrate mix for Aurora is 500μM Kemptide peptide (LRRASLG (SEQ ID NO: 1), Upstate Discovery) in 1ml of water with 35 μCi γ^{33} P-ATP. The substrate mix for GSK3-β is 12.5 μM phospho-glycogen synthase peptide-2 (Upstate Discovery) in 1ml of water with 35 μCi γ^{33} P-ATP. Enzyme and substrate are added to 96 well plates along with 5 μl of various dilutions of the test compound in DMSO (up to 2.5%). The reaction is allowed to proceed for 30 minutes (Aurora) or 3 hours (GSK3-β) before being stopped with an excess of ortho-phosphoric acid (5 μl at 2%). The filtration procedure is as for Activated CDK2/CyclinA assay above.

Please delete the paragraph bridging pages 158-159, and replace it with the following paragraph:

EXAMPLE 50

Measurement of inhibitory activity against Glycogen Synthase Kinase-3 (GSK-3)

GSK3β (human) is diluted to a 10x working stock in 50mM Tris pH 7.5, 0.1mM EGTA, 0.1mM sodium vanadate, 0.1% β-mercaptoethanol, 1mg/ml BSA. One unit equals the incorporation of 1nmol of phosphate per minute phospho-glycogen synthase peptide 2 per minute.

In a final reaction volume of 25μl, GSK3β (5-10 mU) is incubated with 8mM MOPS 7.0, 0.2mM EDTA, 20μM YRRAAVPPSPSLSRHSSPHQS(p)EDEEE (phospho GS2 peptide)

<u>SEQ ID NO: 2)</u>, 10mM MgAcetate and $[\gamma^{-33}P\text{-ATP}]$ (specific activity approx 500cpm/pmol, concentration as required). The reaction is initiated by the addition of Mg²+[$\gamma^{-33}P\text{-ATP}$]. After incubation for 40 minutes at room temperature the reaction is stopped by the addition of 5µl of a 3% phosphoric acid solution. 10µl of the reaction is spotted onto a P30 filter mat and washed 3 times for 5 minutes in 50mM phosphoric acid and once in methanol prior to drying and counting.